Maintenance & Tips
bontempi

Surgical Instruments Care and Maintenance

Rinse and Clean / Sharpen / Tips / Service & Repair

User Guide
bontempi

Instructions and Methods of Use

Cleaning & Sterilization

Surgical Instruments Care and Maintenance

Rinse and Clean / Sharpen / Tips / Service & Repair

Instructions and Methods of Use Cleaning & Sterilization

Bontempi, founded in 1974, is involved in the production and marketing of surgical instruments with the objective of designing instruments to enhance the technical abilities of the individual doctors and professionals towards becoming more efficient and attain better results.

Bontempi's complete product portfolio, of over 10,000 instruments, is positioned as one of the few companies in the world that caters to the **Dental, Plastic, Veterinary, and General surgery** markets.

Globally, Bontempi has three facilities located in Germany, Canada, and Italy, and a number of growing Distributors and Representative offices located all over the world, from Europe to the United States, from the Balkans to the Middle East, Eastern Europe and Asia.

Bontempi instruments are **made in Germany** where the world's best artisans and martensitic stainless steels (AISI 421, AISI 440, AISI 440C) of highest quality can be found.

All products have a **5 year warranty** against manufacturer defects.

Choose quality, choose Bontempi.

Surgical Instruments Care and Maintenance

Rinse and Clean


WHAT? Cleanlact is a concentrated anti-corrosive and bacteriostatic detergent.

WHY? It has an anti-corrosive, bacteriostatic and lubricating action on all surgical instruments; it guarantees perfect maintenance by eliminating eventual stains, halos and browning due to sterilization, disinfection and washing with tap water.

WHEN? Use periodically to extend instruments life. Cleanlact contains a preservative system that prevents bacterial growth.

Cleanlact
DB-10M / 250 cc
DB-50M / 500 cc
DB-1L / 1000 cc


Brush DB-11
155 mm


WHAT? The brushes DB-11 and DB-12 are specially designed for cleaning surgical instruments.

WHY? To remove eventual residue, organic or not, from all surgical instruments, in particular in angled, hinged and knurled parts, without damaging instrument surface. They are autoclavable.

WHEN? Every time, for regular cleaning of all instrumentation and in the event that an instrument is exposed for an extended period to air before being immersed in the cleaning solution.

Brush DB-12
175 mm

Special Rubber for metal
DB-07 Fine
DB-08 Coarse
50 x 40 x 20 mm


WHAT? The Special Rubber for metal is an abrasive eraser specially designed for surgical instruments.

WHY? It is used to remove stains, halos or persistent burnishing from all surgical instruments without damaging the instrument surface.

WHEN? It is necessary when routine cleaning does not remove all stains.

Sharpen

Sharpening Stones

Arkansas N°8 DB-96
100 x 50 x 13 mm

WHAT? Arkansas N° 8 is a natural stone extracted from the mountains of Arkansas: its shape is flat and rectangular.

WHY? It is used to sharpen all cutting parts of surgical instruments.

WHEN? It is necessary to sharpen cutting parts of surgical instruments (scissors, osteotomes, curettes, etc.).

Arkansas Conica N° 299 DB-299
90 mm / Ø 8 mm

WHAT? Arkansas Conica N° 299 is a natural stone extracted from the mountains of Arkansas; it is cone-shaped and cylindrical.

WHY? It is used to manually sharpen internal parts of curettes, scalers, gouges, rongeurs, etc.

WHEN? It is necessary to sharpen the internal parts of surgical instruments in areas that otherwise cannot be reached.

Arkansas N°4 DB-4
100 x 12 x 25 mm

WHAT? Arkansas N° 4 is a natural stone extracted from the mountains of Arkansas; its shape is flat and rectangular.

WHY? It is used to sharpen all cutting parts of surgical instruments.

WHEN? It is used in routine sharpening and finishing of the edge and toe of an instrument.

India N°6 (medium grit) DB-3
Arkansas N°6A (fine texture): DB-6A
100 x 45 mm

WHAT? India stones are quarried from the finest Ozark novaculite (silicon quartz) deposits. Arkansas stones are from the mountains of Arkansas.

WHY? To manually sharpen internal parts of curettes and scalers.

WHEN? India: for sharpening of excessively dull instruments or those requiring re-contouring. Arkansas: for routine sharpening and finishing.

Scharpenoil DB-OL
30 cc

WHAT? Scharpenoil is a non-coloured oil, which is an indispensable complement for manual sharpening of surgical instruments.

WHY? To activate the abrasive properties of all stones (Arkansas, India, etc.), and to reduce the possibility of overheating the instruments during sharpening.

WHEN? Must ALWAYS be used during all sharpening operations.

Lubrication & Testing

WHAT? Rodhoil is a high density lubricating oil of slightly red tint.

WHY? To ensure perfect function and fluid movement of all surgical instruments.

WHEN? It should be applied on all hinges or all other parts which are subject to movement or friction. Its use facilitates elimination of any residue present within the instrument hinges.

Lubricant Rodhoil DB-OLR
50 cc


Teflon Stick for testing DB-15
100 mm / Ø 8 mm

WHAT? Hard plastic test stick that can be used to determine the sharpness of an instrument.

WHY? Use the Bontempi Teflon Stick to test the sharpness of an instrument.


WHEN? If the blade of the instrument runs smoothly over the plastic stick, then it is blunt. A sharp instrument will grab into the stick and removes small fragments of the plastic.

Sterilization

Sterilization Containers

CT-555N / CT-556N
300 x 140 x 40 mm
CT-560N / CT-564N
300 x 140 x 70 mm
CT-565N / CT-569N
300 x 140 x 100 mm

Stainless Steel Cassettes with silicone racks

BS-60 (10 instruments)
185 x 135 x 35 mm
BS-61 (20 instruments)
290 x 185 x 35 mm

Surgical Instruments Care and Maintenance

Rinse & Clean

Immediately after surgery, rinse instruments under warm (not hot) running water. Use only neutral PH detergents.

Ultrasound Cleaning - DO NOT USE for sharp and delicate instruments!

It is recommended for regular instruments such as hand instruments and forceps, following this procedure:

- 1- Sort instruments carefully so as to include only instruments compatible with ultrasonic cleaning.
- 2- Place instruments in open position.
- 3- Make certain sharp edges are not touching other instruments.
- 4- Do not combine different metals (stainless, copper, chrome plated, etc.)
- 5- Change solution frequently to avoid accumulation of microorganisms.
- 6- After rinsing and before sterilization, inspect and dry the instruments thoroughly.

Sharpen

WHY? Surgical Instruments should be kept identical to their original design. Dental and surgical procedures are most effective when using sharp instruments as they reduce hand and wrist fatigue, improve calculus removal, save time, improve tactile sensitivity, and minimize patient discomfort.

WHEN? Instruments should be sharpened lightly after each use; there are two ways to evaluate whether the cutting edge is dull and requires sharpening:

- 1- Visual. The cutting edge should be inspected regularly in a good light (and if possible, under magnification). If the cutting edge is blunt it will be rounded and reflect the light. A dull, non-reflective line indicates sharpness.
- 2- Test Stick. If the blade of the instrument runs smoothly over the testing stick (DB-15), then it is blunt. A sharp instrument will grab into the stick and removes small fragments of the plastic.

HOW? Generally, the instruments should be sharpened following this procedure:

- 1- Place one drop of **Bontempi sharpening oil** (DB-OL) on the **Arkansas stone**. Lubrication improves the movement of the instrument blade over the stone; also, it prevents the metal particles from clogging the stone.
- 2- Hold the instrument in one hand, while applying the stone to the lateral surface angled with the face of the blade.
- 3- Position the stone to contact the heel of the blade and work toward the tip, keeping the stone in contact with the blade throughout the sharpening procedure.
- 4- Move the stone up and down with short strokes, placing more pressure on the down stroke. (Do not move the instrument, keep the instrument still)
- 5- Finish sharpening the instrument with a down stroke; this will prevent a rough edge from forming.
- 6- Evaluate the sharpness with the testing stick (Bontempi code DB-15). If the blade is still dull, re-evaluate the angle of the stone and repeat the sharpening procedure (steps 2-6)

Scalers and Curettes

Sickle Scaler Toe End: the sickle scaler has a pointed tip and, therefore, the stone is held straight as it nears the tip.


Curette Toe End: the curette has a rounded toe, so the position of the stone is adapted around the rounded cross-section. With both types of instrument, always finish on a down stroke to remove any flash of metal.

Care of Sharpening Stones

After use, wipe the stone with a clean cloth to remove metal particles. Then, clean the stone by scrubbing or using ultrasound to remove lubricant before sterilization. After sterilization, lubricate before each use. Be sure to use entire stone to prevent "grooving".

Service & Repair

- **Sharpening**
(Scissor, Bone Curette, Bone Cutter, Elevator, Rasp, Rongeur, Osteotome, Scalers, Periosteals, Knives)
- **Tungsten-Carbide Inserts Replacement**
(Needle Holders)
- **Alignment**
(Hemostat, Forceps)
- **Polishing**
(All instruments)
- **Personalized Etching**
(All instruments)
- **Repairs**
(Instrument to be assessed)

Lubrication:
Regular use of the lubricant "Rodhoil" (DB-OLR) and "Cleanlact" (DB-50M), will prevent rust, corrosion, and stiff joints.

Ultrasound Cleaning:
Is not recommended for scissors, needle holders, bone forceps and all tungsten carbide (TC) inserts.

Sharpening:
Regular sharpening will enhance the life of your cutter and will improve the quality of your work.

Sharp & Delicate Instruments:
Scissors, needle holders or bone forceps should be cleaned manually by using Bontempi cleaning brushes (DB-11 and DB-12).

Tungsten carbide Instruments:
TC instruments are more sensitive to chemicals and require special care. they should never be exposed to chemical substances or to any other corrosive chemicals.

Instructions and Methods of Use

Cleaning and Sterilization

All Bontempi instruments are hand crafted in Germany from grade 440 Rockwell surgical steel, in order to exceed surgeon specifications for ergonomics and duration. Stainless steel is mainly used for the production of all instruments that do not require sharpening, while carbon steel, a special alloy, retains characteristics that are more suitable for cut and less for corrosion. Due to the high quality of our surgical products, Bontempi is proud to offer a five-year guarantee covering all manufacturing defects. In order to ensure proper use of our instruments, we suggest you read the following instructions on product care to ensure their maximum efficiency and duration.

FIRST USE

The instrument is supplied in a "NON STERILE" and imbued with a special liquid oily in cycles, parts and pivoting in those screwed. Therefore, before use, you must good wash the instruments using neutral detergent and / or degreasers suitable for steel, so as not to cause unpleasant complications with the instrument (red spots, induration pivoting shares, etc.).

CLEANING

Before being sent to destination, all our instruments are inspected, lubricated, and individually packed. It is necessary to wash and cleanse each instrument carefully prior to each usage. After having used the instruments, brush off any excess particles using nylon bristles (never use steel bristles) and rinse them individually. In the event of a potentially infected instrument, soak the instrument in a disinfectant detergent solution for at least 10 minutes. Cleaning in an ultrasound basin is not recommended for instruments with cutting edges (scissors, bistouries, bone forceps, scalpels, etc.) and instruments with Tungsten carbide insertions (TC): they may chip, break or corrode. Cutting instruments of different materials (chromate, copper, aluminium, and stainless steel) should always be cleaned and stored separately. After the cleaning, make sure that all the instruments are perfectly dry. We recommend re-lubricating the instruments with surgically approved products only. Never use industrial oils or lubricants. After having cleaned, rinsed and lubricated the instruments the sterilizing phase can start.

STERILIZATION

Sterilization is a procedure which removes all pathogenic and non pathogenic micro-organisms; it should be performed with outmost care. Remember that instruments sterilization cannot substitute cleaning. These are the most common and effective sterilization techniques:

Autoclave (Steam Sterilization)

Instruments are treated with high pressure water steam.

Recommended standard time and temperature*: 20 minutes, at a temperature of 134°C (286° F).

Sterilizing agent: distilled water.

Chemiclave

Instruments are treated with high pressure, with the help of chemical steams. Instruments must be dried before sterilization.

Recommended standard time and temperature*: 20 minutes, at a temperature of 130°C (266° F).

Sterilizing agent: special chemical solutions.

Dry Sterilization

Instruments sterilization takes place by means of hot air. Instruments must be dried before sterilization.

Recommended standard time and temperature*: from 60 to 90 minutes, at a temperature of 150°/170° C (300°/320° F).

At present, this sterilization technique is not recommended, as it doesn't guarantee reliable sterilization and gradually causes damages such as surface blackening, loss of polish and, especially for cutting instruments, dullness.

* Recommended times and temperatures may vary according to the nature and loading capacity of your equipment and do not take heating times into consideration. Always comply with the manufacturer's instructions. The above mentioned recommendations concerning maintenance, cleaning and sterilization of instruments aim at the safeguard of both patient and operator. Remember that if these procedures are carried out correctly, this will enhance the Quality of Bontempi instruments, as well as prolong their effectiveness over time.

WARNINGS

1-Used, damaged and oxidized tools should not be used because they are no longer able to perform their function. Please note that some types of damages (corrosion, rust and spots) are transmitted to the instruments intact.

2-Avoid putting oxidized or rusty instruments in sterilizing or disinfectant solutions, as other instruments could be attacked.

3-Bontempi's product warranty against manufacturer defects automatically expires in the cases of improper care, maintenance, and use of the products.

4-The tools provided by Bontempi are subject to Directive 93/42/EEC and therefore be discarded by the regional force.

Each Bontempi instrument is produced and controlled with the utmost care and destined exclusively for professionals and for the use which they have been devised and realized.

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